

LIFEPAK CR2 Defibrillator Features Guide

[AED Machine](#) | [LIFEPAK CR2 Defibrillator](#)

At AED Advantage we supply the LIFEPAK CR2—an innovative yet easy-to-use AED trusted by professionals and first responders. Its advanced technology (including CPRInsight™), clear instructions and Wi-Fi monitoring make life-saving interventions as simple as possible.

With user-friendly design, durable construction and child mode, the CR2 meets the needs of workplaces, schools and public spaces. Always connected via **Wi-Fi** for real-time monitoring, it ensures your AED is ready when you need it most.



LIFEPAK CR2 Defibrillator Features Guide

[AED Machine](#) | [LIFEPAK CR2 Defibrillator](#)



CPRInsight™ Technology

The **CPRInsight™** feature allows chest compressions to continue uninterrupted during rhythm analysis, increasing hands-on time and improving survival rates. It automatically differentiates between patient movement and the ECG signal.

- Analyzes heart rhythm during compressions
- Minimizes pauses for higher CPR quality
- Available on Wi-Fi enabled models



User-Friendly Interface

Designed for lay responders and professionals alike, the CR2 guides users through each step with clear voice prompts and intuitive graphics.

- Calm voice prompts reduce stress
- Simple icons guide pad placement
- Helps untrained bystanders act confidently



LIFEPAK CR2 Defibrillator Features Guide

[AED Machine](#) | [LIFEPAK CR2 Defibrillator](#)



Wi-Fi Connectivity

Integrated **Wi-Fi** links the AED to AED Advantage's remote monitoring platform so your CR2 stays ready without manual inspections.

- Real-time status checks
- Alerts for battery and pad expiry
- Ideal for managing multiple units



Durable & Portable Design

Despite its compact size, the CR2 is rugged enough for demanding environments and easy to transport between locations.

- Sturdy housing withstands bumps
- Lightweight with carrying handle
- Long-life battery



LIFEPAK CR2 Defibrillator Features Guide

[AED Machine](#) | [LIFEPAK CR2 Defibrillator](#)



Child Mode

A dedicated **child** mode allows rescuers to switch to paediatric energy levels without changing electrodes, simplifying mixed-age environments.

- Instantly adjusts energy for children
- No need for separate pads
- Suitable for schools and public spaces